

Product datasheet for **AP05381HR-N**

Rat IgG2c Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2c Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1:10000. Western Blot: 1:1000 - 1:10000.
Reactivity:	Rat
Host:	Goat
Isotype:	IgG
Formulation:	PBS containing 0,1% ProclinTM300 and 0,2% BSA Label: HRP State: Liquid purified IgG Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)