

## Product datasheet for **AP05382HR-N**

### Rat IgM Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1:10,000 - 1:100,000. Western Blot: 1:1000 - 1:10,000.
Reactivity:	Rat
Host:	Goat
Isotype:	IgG
Formulation:	Phosphate buffered saline containing 0.1% Proclin™ 300, 0.2% Bovine Serum Albumin Label: HRP State: Liquid purified IgG Label: Horseradish Peroxidase
Concentration:	lot specific
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)