

## Product datasheet for **AP05382PU-N**

### Rat IgM Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM Goat Polyclonal Antibody
Applications:	ELISA, FC, WB
Recommended Dilution:	ELISA 1ug/ml - 10ug/ml, Western Blotting, Flow Cytometry
Reactivity:	Rat
Host:	Goat
Isotype:	IgG
Formulation:	Phosphate buffered saline pH7.2 0.09% Sodium Azide, State: Purified State: Purified IgG - liquid
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. 18 months from date of despatch.



[View online »](#)