

## Product datasheet for **AP05383AP-N**

### Mouse IgG (adsorbed) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (adsorbed) Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA:</b> 1/200-1/2500. <b>Western Blot:</b> 1/100-1/500.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Whole Mouse IgG
Isotype:	IgG
Formulation:	50mM HEPES pH 7.1, 0.1M NaCl, 0.1M MgCl <sub>2</sub> , 0.1M ZnCl <sub>2</sub> Label: AP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)