

## Product datasheet for **AP05383FC-N**

### Mouse IgG (adsorbed) Goat Polyclonal Antibody


#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (adsorbed) Goat Polyclonal Antibody
Applications:	FC
Recommended Dilution:	<b>Flow Cytometry:</b> Use 50 µl of 1/50-1/200 diluted antibody to label 1x10 <sup>6</sup> cells in 100 µl.
Reactivity:	Mouse
Host:	Goat
Isotype:	IgG
Formulation:	PBS pH 7.2 Label: FITC State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)

## Product images:



Rabbit peripheral blood granulocytes stained with Mouse anti Rabbit CD44 (Cat.-No SM607) followed by Goat anti Mouse IgG:FITC (Cat.-No AP05383FC-N)