

Product datasheet for **AP05383HR-N**

Mouse IgG (adsorbed) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (adsorbed) Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/10000-1/100000. Western Blot: 1/1000-1/10000.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Antisera to Mouse IgG were raised by repeated immunisations of Goats with highly purified antigen.
Isotype:	IgG
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.1% Proclin TM 300 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)