

Product datasheet for **AP05384AP-N**

Mouse IgG1 Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG1 Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/2000-1/4000. Western Blot. Immunohistochemistry on Frozen Sections. Immunohistochemistry on Paraffin Sections.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG1 paraproteins.
Isotype:	IgG
Formulation:	50mM TRIS, 1mM MgCl ₂ , pH 8.0 Label: AP State: Liquid purified IgG fraction from whole serum Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)