

Product datasheet for **AP05384FC-N**

Mouse IgG1 Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG1 Goat Polyclonal Antibody
Applications:	FC, IF, IHC
Recommended Dilution:	Flow Cytometry: Use 10 µl of neat-1/10 diluted antibody to label 1 ⁰ cells in 100 µl. Immunohistochemistry on Frozen Sections. Immunocytochemistry: neat-1/10.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG1 paraproteins.
Isotype:	IgG
Formulation:	PBS Label: FITC State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use. Avoid repeated freezing and thawing.



[View online »](#)