

Product datasheet for **AP05384HR-N**

Mouse IgG1 Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG1 Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/4000-1/8000. Western Blot. Immunohistochemistry on Frozen Sections. Immunohistochemistry on Paraffin Sections.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG1 paraproteins.
Isotype:	IgG
Formulation:	PBS, pH 7.4 Label: HRP State: Liquid purified IgG fraction Stabilizer: 50% Glycerol Preservative: 0.01% Thiomersal Label: Horseradish Peroxidase
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)