

Product datasheet for AP05384RP-N

Mouse IgG1 Goat Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Mouse IgG1 Goat Polyclonal Antibody

Applications: FC

Recommended Dilution: Flow Cytometry: 1/20 - 1/30.

Reactivity: Mouse **Host:** Goat

Immunogen: Mouse IgG1 paraproteins

Isotype: IgG

Formulation: Phosphate buffered saline pH7.4 containing 0,09% sodium azide and sucrose

Label: PE

State: Liquid purified IgG Label: R.Phycoerythrin (RPE)

Concentration: lot specific

Conjugation: PE

Storage: Store the antibody undiluted at 2-8°C.

DO NOT FREEZE!

This product is photosensitive and should be protected from light.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com