

Product datasheet for **AP05385AP-N**

Mouse IgG2a Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/2000 - 1/4000.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2a paraproteins
Isotype:	IgG
Formulation:	50mM TRIS, 1mM MgCl ₂ , pH8,0 containing 0,09% sodium azide and 50% glycerol State: Liquid purified IgG Label: Alkaline Phosphatase
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)