

## Product datasheet for **AP05385FC-N**

### Mouse IgG2a Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a Goat Polyclonal Antibody
Applications:	FC, IF, IHC
Recommended Dilution:	Flow Cytometry: 1/10. Immunofluorescence: 1/10.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2a paraproteins
Isotype:	IgG
Formulation:	PBS pH7.4 containing 0,09% sodium azide Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)