

Product datasheet for **AP05385RP-N**

Mouse IgG2a Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a Goat Polyclonal Antibody
Applications:	FC
Recommended Dilution:	Flow Cytometry: 1/20 - 1/30.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2a paraproteins
Isotype:	IgG
Formulation:	PBS, pH 7,4 containing 0,09% sodium azide and Sucrose Label: PE State: Liquid purified IgG Label: R.Phycoerythrin (RPE)
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light.



[View online »](#)