

## Product datasheet for **AP05386AP-N**

### Mouse IgG2b Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2b Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA:</b> 1/2000-1/4000. <b>Western Blot.</b> <b>Immunohistochemistry on Frozen Sections.</b> <b>Immunohistochemistry on Paraffin Sections.</b>
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2b paraproteins.
Isotype:	IgG
Formulation:	TRIS buffered saline
	Label: AP State: Liquid purified IgG fraction Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)