

## Product datasheet for **AP05386RP-N**

### Mouse IgG2b Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2b Goat Polyclonal Antibody
Applications:	FC
Recommended Dilution:	<b>Flow Cytometry:</b> Use 10 µl of 1/20-1/30 diluted antibody to label 1 <sup>0</sup> cells in 100 µl.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG2b paraproteins.
Isotype:	IgG
Formulation:	PBS, pH 7.4 containing 0.09% Sodium Azide and Sucrose Label: PE State: Liquid purified IgG fraction Label: R. Phycoerythrin (RPE)
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light.



[View online »](#)