

Product datasheet for **AP05387AP-N**

Mouse IgG2c Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2c Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000 - 1/10000.
Reactivity:	Mouse
Host:	Goat
Isotype:	IgG
Formulation:	50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂ containing 0.09% Sodium Azide, 0.2% Bovine Serum Albumin State: Liquid purified IgG Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)