

## Product datasheet for **AP05389HR-N**

### Mouse IgA Goat Polyclonal Antibody

#### Product data:

|                       |   |
|-----------------------|---|
| Product Type:         | Secondary Antibodies  |
| Product Name:         | Mouse IgA Goat Polyclonal Antibody  |
| Applications:         | ELISA, IHC, WB  |
| Recommended Dilution: | Immunohistochemistry on frozen and paraffin sections.<br>ELISA: 1/4000 - 1/8000.<br>Western Blot.   |
| Reactivity:           | Mouse   |
| Host:                 | Goat  |
| Immunogen:            | Mouse IgA paraproteins  |
| Isotype:              | IgG   |
| Formulation:          | Phosphate buffered saline pH7.4 containing 0.01% Thiomersal and 50% Glycerol<br>Label: HRP<br>State: Liquid purified Ig fraction<br>Label: Horseradish Peroxidase |
| Purification:         | Affinity chromatography   |
| Conjugation:          | HRP   |
| Storage:              | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.                                 |



[View online »](#)