

## Product datasheet for **AP05389PU-N**

### Mouse IgA Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA Rabbit Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Reactivity:	Mouse, Rat
Host:	Rabbit
Immunogen:	Mouse IgA paraproteins
Isotype:	IgG
Formulation:	Borate buffered saline pH8.2 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)