

Product datasheet for **AP05390AP-N**

Mouse IgM Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgM Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry on frozen and paraffin sections. ELISA: 1/2000 - 1/4000. Western Blot.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgM paraproteins
Isotype:	IgG
Formulation:	50mM TRIS, 1mM MgCl ₂ , pH8.0 containing 0.09% Sodium Azide (NaN ₃) and 50% Glycerol State: Liquid purified IgG Label: Alkaline phosphatase
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)