

## Product datasheet for **AP05391FC-N**

### Human IgD (heavy chain) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD (heavy chain) Goat Polyclonal Antibody
Applications:	FC, IF, IHC
Recommended Dilution:	<b>Flow Cytometry:</b> Use 10 µl of Neat-1/10 diluted antibody to label 10 <sup>6</sup> cells in 100 µl. <b>Immunocytochemistry:</b> Neat-1/10. <b>Immunohistochemistry on Frozen Sections.</b>
Reactivity:	Human
Host:	Goat
Immunogen:	Human IgD paraproteins.
Isotype:	IgG
Formulation:	PBS Label: FITC State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)