

## Product datasheet for **AP05392PU-N**

### Human IgD (F(ab)2 Fragment) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD (F(ab)2 Fragment) Rabbit Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA
Reactivity:	Human
Host:	Rabbit
Immunogen:	Human IgD
Isotype:	IgG
Formulation:	Borate buffered saline pH 8.2 containing 0.09% Sodium Azide State: Purified State: Liquid F(ab)2 fragments
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)