

## Product datasheet for **AP05393AP-N**

### Human IgE Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/1000-1/10000. Western Blot: 1/1000-1/5000.
Reactivity:	Human, Monkey
Host:	Goat
Isotype:	IgG
Formulation:	50mM HEPES containing 0.1M NaCl, 1mM MgCl <sub>2</sub> , 0.1mM ZnCl <sub>2</sub> with 0.09% Sodium Azide as preservative and 0.2% BSA as stabilizer. Label: AP State: Liquid purified IgG fraction. Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)