

## Product datasheet for AP05393BT-N

## Human IgE Goat Polyclonal Antibody

## **Product data:**

Product Type:	Secondary Antibodies
Product Name:	Human IgE Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/10,000-1/100,000. Western Blot: 1/10,000-1/100,000.
Reactivity:	Human, Monkey
Host:	Goat
Immunogen:	Native human IgE
lsotype:	IgG
Formulation:	PBS Label: Biotin State: Liquid purified IgG fraction. Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US