

Product datasheet for **AP05393FC-N**

Human IgE Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Goat Polyclonal Antibody
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	Human, Monkey
Host:	Goat
Immunogen:	Native human IgE
Isotype:	IgG
Formulation:	PBS Label: FITC State: Liquid purified IgG fraction. Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light.



[View online »](#)