

Product datasheet for **AP05393HR-N**

Human IgE Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10000-1/100000.
Reactivity:	Human, Monkey
Host:	Goat
Immunogen:	Native Human IgE
Isotype:	IgG
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.1% Pro-Clean™400 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Note:	Antisera to human IgE were raised by repeated immunisations of goats with highly purified antigen



[View online »](#)