

Product datasheet for **AP05394PU-N**

Rabbit IgG Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG Sheep Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/1000. Immunohistochemistry on Frozen Sections.
Reactivity:	Rabbit
Host:	Sheep
Immunogen:	Purified Rabbit IgG
Isotype:	IgG
Formulation:	Glycine Buffered Saline pH 7.4 containing 0.09% Sodium Azide, 0.1% E-amino-n-Caproic Acid, 0.01% Benzamidine, 1mM EDTA State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)