

Product datasheet for **AP05395AP-N**

Sheep/Goat IgG Donkey Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Sheep/Goat IgG Donkey Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	Western blot. ELISA: 1/2000-1/5000. Immunohistochemistry on Frozen Sections.
Reactivity:	Goat, Sheep
Host:	Donkey
Immunogen:	Purified IgG from a Sheep serum pool
Isotype:	IgG
Formulation:	Tris Buffered Saline, pH 8.0 containing 0.01% Sodium Azide, 1% BSA, 50% Glycerol. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! Should this product contain a precipitate we recommend microcentrifugation before use.



[View online »](#)