

Product datasheet for **AP05395FC-N**

Sheep/Goat IgG Donkey Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Sheep/Goat IgG Donkey Polyclonal Antibody
Applications:	IHC
Recommended Dilution:	Immunohistochemistry on Frozen Sections: 1/100-1/200.
Reactivity:	Goat, Sheep
Host:	Donkey
Immunogen:	Sheep IgG
Isotype:	IgG
Formulation:	PBS containing 0.09% Sodium Azide and 1% BSA. Label: FITC State: Liquid purified IgG fraction Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)