

Product datasheet for **AP05395HR-N**

Sheep/Goat IgG Donkey Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Sheep/Goat IgG Donkey Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/12800-1/25600. Western Blot: 1/20000. Immunohistochemistry on Paraffin Sections. Immunohistochemistry on Frozen Sections: 1/50-1/100.
Reactivity:	Goat, Sheep
Host:	Donkey
Immunogen:	Sheep IgG
Isotype:	IgG
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction. Stabilizer: 1% BSA Preservative: 0.05% Proclin™ 300, 0.002% Methylisothiazolone and 0.002% Bromonitrodioxane Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)