

Product datasheet for **AP05413PU-N**

E. coli K88AB Pili Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	ELISA.
Reactivity:	Escherichia coli
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E. coli strain 08:K87:K88ab was used.
Specificity:	This antibody reacts to <i>Escherichia coli</i> K88B Pili
Formulation:	Phosphate buffered saline containing 0.09% Sodium Azide (NaN ₃) State: Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Ammonium sulphate precipitation
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Escherichia coli, EHEC, EPEC, ETEC



[View online »](#)