

## Product datasheet for **AP05440SU-N**

### Growth Hormone (GH1) Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	Radioimmunoassays. ELISA.
Reactivity:	Human
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Purified human growth hormone.
Specificity:	This antibody is specific for human growth hormone.
Formulation:	State: Serum State: Liquid serum containing 0.09% Sodium Azide (NaN <sub>3</sub> )
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	growth hormone 1
Database Link:	<a href="#">Entrez Gene 2688 Human P01241</a>
Background:	Growth hormone is a small peptide hormone released by the anterior pituitary under the control of the hypothalamus. It acts to control body growth and metabolic activity.
Synonyms:	Pituitary growth hormone, Growth hormone 1, HGH



[View online »](#)