

Product datasheet for **AP05448PU-N**

Human Kappa Light Chain Sheep Polyclonal Antibody

Product data:

| | |
|-----------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, IHC |
| Recommended Dilution: | Immunohistochemistry on frozen sections: 1/50 - 1/200. ELISA. |
| Reactivity: | Human |
| Host: | Sheep |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | Native |
| Specificity: | This antibody reacts to kappa light chain. |
| Formulation: | Phosphate buffered saline containing 0.09% Sodium Azide (NaN ₃) State: Purified State: Liquid Ig fraction |
| Concentration: | lot specific |
| Purification: | Ammonium sulphate precipitation |
| Conjugation: | Unconjugated |
| Storage: | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |



[View online »](#)