

Product datasheet for **AP05497SU-N**

Progesterone Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Broad
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	HS-HMW conjugate
Specificity:	Exhibits some cross reactivity with hydroxycortisone, but minimal reactivity with 17-a-hydroxyprogesterone, 5-a-pregnane-3,20-dione, corticosterone, 17-a-hydroxypregnenolone, 5-a-dihydroxyprogesterone.
Formulation:	State: Serum State: Liquid serum containing 0.09% Sodium Azide
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	20-dione, pregn-4-ene-3



[View online »](#)