

Product datasheet for **AP05508SU-N**

Propranolol Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Propranolol-hemisuccinate conjugated to BSA.
Specificity:	This antibody recognises Propranolol, a non-selective beta blocker used for the treatment of hypertension and Angina pectoris.
Formulation:	State: Serum State: Liquid serum
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)