

Product datasheet for **AP05548HR-N**

Goat IgA Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat IgA Rabbit Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000-1/100,000.
Reactivity:	Goat
Host:	Rabbit
Immunogen:	Purified goat IgA.
Isotype:	IgG
Formulation:	Phosphate buffered saline containing 0.2% BSA as stabilizer and 0.1% Pro-clean(TM)400 as preservative. Label: HRP State: Liquid purified IgG fraction. Label: Horseradish Peroxidase
Concentration:	lot specific
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)