

Product datasheet for **AP05549HR-N**

Goat IgM Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat IgM Rabbit Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000-1/100,000.
Reactivity:	Goat
Host:	Rabbit
Immunogen:	Purified Goat IgM
Isotype:	IgG
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.1% Pro-clean™400 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)