

Product datasheet for **AP05549PU-N**

Goat IgM Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat IgM Rabbit Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA: 1/100 - 1/1000. Immunodiffusion.
Reactivity:	Goat
Host:	Rabbit
Immunogen:	Purified goat IgM
Isotype:	IgG
Formulation:	Phosphate buffered saline containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)