

## Product datasheet for **AP05550AP-M**

### Rabbit IgA Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgA Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/100-1/1000.
Reactivity:	Rabbit
Host:	Goat
Immunogen:	Purified Rabbit IgA
Isotype:	IgG
Formulation:	50mM HEPES, 0.1M NaCl, 1mM MgCl <sub>2</sub> , 0.1mM ZnCl <sub>2</sub>
	Label: AP
	State: Liquid purified Ig fraction
	Stabilizer: 0.2% BSA
	Preservative: 0.09% Sodium Azide
	Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>
Note:	<b>Antiserum Preparation:</b> Antisera to Rabbit IgA were raised by repeated immunisation of goats with highly purified antigen.



[View online »](#)