

## Product datasheet for **AP05550FC-N**

### Rabbit IgA Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgA Goat Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	Immunohistochemistry on frozen sections: 1/200 - 1/2000. Flow Cytometry.
Reactivity:	Rabbit
Host:	Goat
Immunogen:	Purified rabbit IgA
Isotype:	IgG
Formulation:	Phosphate buffered saline containing 0.09% Sodium Azide, 0.2% Bovine Serum Albumin Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)