

Product datasheet for **AP05550HR-N**

Rabbit IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgA Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10000-1/100000.
Reactivity:	Rabbit
Host:	Goat
Immunogen:	Purified rabbit IgA
Isotype:	IgG
Formulation:	PBS containing 0.1% Pro-clean(TM) 400, 0.2% BSA as stabilizer Label: HRP State: Liquid purified IgG fraction Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)