

Product datasheet for **AP05550PU-N**

Rabbit IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgA Goat Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA: 1/100 - 1/1000. Immunodiffusion.
Reactivity:	Rabbit
Host:	Goat
Immunogen:	Purified rabbit IgA
Isotype:	IgG
Formulation:	Phosphate buffered saline pH7.2 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)