

## Product datasheet for **AP05966SU-N**

### Aldosterone Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Reactivity:	Human
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Specificity:	This antibody is specific for Human Aldosterone.
Formulation:	State: Serum State: Liquid Serum Preservative: 0.09% Sodium Azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Aldosterone is a steroid hormone which is produced by the adrenal gland. It has a number of functions relating to the kidney nephron including the balance of sodium and potassium, and increasing blood pressure. It is involved in two disorders: Addisons disease, where it is reduced, and Conn Syndrome, where it is found in increased amounts.



[View online »](#)