

## Product datasheet for **AP05970SU-N**

### Paraquat Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Paraquat N-valeric acid derivative (N- (4-carboxybutyl)-N'-methyl-bipyridilium) conjugated to keyhole limpet haemocyanin (KLH).
Specificity:	This antibody recognises Paraquat, a bipyridyl herbicide. It is widely used but is extremely toxic to humans. The chemical name is N,N'-Dimethyl-4,4'-bipyridinium dichloride.
Formulation:	State: Serum State: Liquid serum
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Methyl viologen



[View online »](#)