

Product datasheet for **AP05973SU-N**

Triiodothyronine (T3) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	Radioimmunoassays.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Triiodothyronine (T3) derivative coupled to Keyhole Limpet Hemocyanon (KLH)
Specificity:	This antibody recognises Triiodothyronine in radioimmunoassays. This antibody has not been cross-adsorbed and may react with related compounds.
Formulation:	State: Serum State: Liquid serum
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)