

Product datasheet for **AP06279PU-M**

PAX5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western blot: 1/500-1/1000.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, corresponding to the N-terminus of Human Pax-5.
Specificity:	This antibody detects endogenous levels of Pax-5 protein. (region surrounding Asp2)
Formulation:	Phosphate buffered saline (PBS), pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction Preservative: 15 mM sodium azide
Concentration:	1.0 mg/ml
Purification:	Affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 40 kDa
Gene Name:	paired box 5
Database Link:	<u>Entrez Gene 5079 Human Q02548</u>




[View online »](#)

Background:

The Pax gene family of nuclear transcription factors is comprised of nine family members that function during embryogenesis to regulate the temporal and position-dependent differentiation of cells. In addition, the family is involved in a variety of signal transduction pathways in the adult organism. Mutations in the Pax family of proteins have been linked to disease and cancer in humans. For example, the human Pax-5 gene encodes a B cell lineage specific protein called B cell specific activator protein or BSAP, which is expressed in pro-B, pre-B and mature B lymphocytes but not in plasma cells. BSAP functions to regulate not only B cell development, but also influences the balance between immunoglobulin secretion and B cell proliferation. Overexpression of BSAP has been implicated in cellular transformation, and in the case of small lymphocytic lymphomas with plasmacytoid differentiation, a t(9;14) (p13;q32) translocation resulting in the deregulation of Pax-5 gene expression has been detected. The gene which encodes Pax-5 maps to human chromosome 9p13.

Synonyms:

Pax-5, BSAP

Product images:


Western blot (WB) analysis of Pax-5 antibody in extracts from A549 cells.