

## Product datasheet for **AP06383PU-S**

### **DARPP32 (PPP1R1B) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	<b>Western blot:</b> 1/500-1/1000. <b>Immunohistochemistry on paraffin sections:</b> 1/50-1/200.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, corresponding to amino acids 35-82 of Human DARPP-32.
Specificity:	This antibody detects endogenous levels of DARPP-32 protein. (region surrounding Pro69)
Formulation:	Phosphate buffered saline (PBS), pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.05% sodium azide
Concentration:	1.0 mg/ml
Purification:	Affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 32 kDa
Gene Name:	protein phosphatase 1 regulatory inhibitor subunit 1B
Database Link:	<a href="#">Entrez Gene 84152 Human</a> <a href="#">Q9UD71</a>




[View online »](#)

**Background:**


Dopaminergic signaling pathways, which are essential for multiple brain functions, are abnormal in several neurological disorders, such as schizophrenia, Parkinson's disease and drug abuse. DARPP-32 (for dopamine and adenosine 3',5'-monophosphate-regulated phosphoprotein of 32 kDa) is abundant in neurons that receive dopaminergic input. Activation of PKA and the consequent phosphorylation of DARPP-32 on threonine occurs in response to dopamine acting upon D1-like receptors. Dopamine interaction with D2-like receptors results in the inhibition of PKA activation, the activation of protein phosphatase 2B and the consequent dephosphorylation of DARPP-32. Neurotransmitters other than dopamine may also be able to stimulate the phosphorylation or dephosphorylation of DARPP-32. Phosphorylated DARPP-32 is a potent inhibitor of PP-1.

**Synonyms:**

PPP1R1B, DARPP-32, FLJ20940

**Product images:**

Western blot (WB) analysis of DARPP-32 antibody in extracts from 293 cells treated with EGF 200ng/ml 30' or HeLa cells treated with EGF 200ng/ml 5'.



Immunohistochemistry (IHC) analyzes of DARPP-32 antibody in paraffin-embedded human brain tissue.