

Product datasheet for **AP06443PU-M**

Nkx3.1 (NKX3-1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Western blot: 1/500-1/1000. Immunohistochemistry on paraffin sections 1/50-1/200.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, corresponding to the N-terminal of Human Nkx-3.1.
Specificity:	This antibody detects endogenous levels of NKX3.1 protein. (region surrounding Leu2)
Formulation:	Phosphate buffered saline (PBS), pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction (> 95% pure by SDS-PAGE). Preservative: 0.05% Sodium Azide
Concentration:	1.0 mg/ml
Purification:	Affinity Chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Predicted Protein Size:	~ 38 kDa
Gene Name:	NK3 homeobox 1
Database Link:	Entrez Gene 4824 Human Q99801




[View online »](#)

Background:

The homeobox gene Nkx-3.1 is the human homolog of *Drosophila* bagpipe, which, in conjunction with tinman, determines cell fate in the dorsal mesoderm. In mammalian species, Nkx-3.1 is predominantly expressed in prostate, and it regulates prostate development in response to sonic hedgehog (Shh) signaling by exerting growth-suppressive and differentiating effects on prostatic epithelium. Nkx-3.1 is also expressed at lower levels in other tissues, including the heart and gut, in a Shh independent manner, where it plays a role in regulating proliferation of glandular epithelium and in the formation of ducts in prostate and minor salivary glands. Nkx-3.1 preferentially binds the TAAGTA sequence, which has not been reported for any other NK class homeoprotein. The human Nkx-3.1 gene is located on chromosome 8q21, which frequently undergoes a loss of heterozygosity, and although Nkx-3.1 is not a tumor suppressor gene, it may be a useful marker for benign and malignant prostate epithelium.

Synonyms:

NKX3.1, NKX3A

Product images:


Immunohistochemistry (IHC) analyzes of NKX3.1 antibody in paraffin-embedded human breast carcinoma tissue.