

Product datasheet for **AP06467PU-N**

TRADD Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Western blot: 1/500-1/1000. Immunohistochemistry on paraffin sections: 1/50-1/200.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody detects endogenous levels of TRADD protein. (region surrounding Leu282)
Formulation:	Phosphate buffered saline (PBS), pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction Preservative: 15 mM sodium azide
Concentration:	1.0 mg/ml
Purification:	Affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 35 kDa
Gene Name:	TNFRSF1A associated via death domain
Database Link:	Entrez Gene 8717 Human Q15628
Background:	TRADD (TNF R1 associated death domain) is a 34kDa protein that associates with the death domains of both FAS and TNFR1 via their DDH regions. Over expression of TRADD leads to NFkB activation and apoptosis in the absence of TNF.
Synonyms:	MGC11078




[View online »](#)

Product images:



Western blot (WB) analysis of TRADD antibody (Cat.-No.: AP06467PU-N) at 1/500 dilution Lane 1:HEK293T whole cell lysate Lane 2:Mouse heart tissue lysate Lane 3:Rat liver tissue lysate



Immunohistochemistry (IHC) analyzes of TRADD antibody (Cat.-No.: AP06467PU-N) in paraffin-embedded human breast carcinoma tissue.