

## Product datasheet for **AP06681PU-N**

### Cleavage stimulation factor 2 (CSTF2) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	<b>Western blot:</b> 1/500-1/1000. <b>Immunohistochemistry on paraffin sections:</b> 1/50-1/200. <b>Immunofluorescence:</b> 1/50-1/200.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, corresponding to amino acids N-terminus of Human CstF-64.
Specificity:	This antibody detects endogenous levels of CstF-64 protein. (region surrounding Arg46)
Formulation:	Phosphate buffered saline (PBS), pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction Preservative: 0.05% sodium azide
Concentration:	1.0 mg/ml
Purification:	Affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 64 kDa
Gene Name:	cleavage stimulation factor subunit 2
Database Link:	<a href="#">Entrez Gene 108062 Mouse</a> <a href="#">Entrez Gene 1478 Human</a> <a href="#">P33240</a>




[View online »](#)

**Background:**

The serine threonine protein kinase encoded by the AKT1 gene is catalytically inactive in serum starved primary and immortalized fibroblasts. AKT1 and the related AKT2 are activated by platelet derived growth factor. The activation is rapid and specific. In the developing nervous system AKT is a critical mediator of growth factor induced neuronal survival. Survival factors can suppress apoptosis in a transcription independent manner by activating the serine/threonine kinase AKT1, which then phosphorylates and inactivates components of the apoptotic machinery. Multiple alternatively spliced transcript variants have been found for this gene (referenced from entrez gene).

**Synonyms:**

CSTF 64 kDa subunit, CstF-64, CSTF64, CF-1 64 kDa subunit

**Product images:**

Western blot (WB) analysis of CstF-64 antibody in extracts from 293 cells.