

Product datasheet for **AP08070PU-N**

Opioid Receptor (pan) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western Blot (1/500-1/1000).
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized non-phosphopeptide derived from Mouse Opioid Receptor around the phosphorylation site of Serine 375 (H-P-SP-T-A).
Specificity:	This antibody detects endogenous levels of total Opioid Receptor protein.
Formulation:	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide and 50% Glycerol. State: Aff - Purified State: Liquid purified Ig fraction.
Purification:	Immunoaffinity Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.



Product images:




Figure 1. Western blot analysis of extract from PC12 cells, using Opioid Receptor antibody AP08070PU (Lane 1 and 2).